

GX Microscopes

product datasheet

XSZ Series

MultiViewing
Biological
Microscopes




Specifications					
		Model:	XSZ N204	XSZ N304	XSZ 510
Eyepieces	Widefield	WF10X (18mm FOV)	STD	STD	STD
		WF16X (11mm FOV)	OPT	OPT	OPT
Objectives	Achromatic	4X/0.10	STD	STD	STD
		10X/0.25	STD	STD	STD
		40X(s)/0.65	STD	STD	STD
		100X(s)/1.25 Oil	STD	STD	STD
Nosepiece	Locating on ball bearing	Quadruple	STD	STD	STD
Additional Viewing Ports	Mounting	On Main Body	2 sets Binocular	1 set Binocular 2 sets Monocular	1 set Binocular
		On Separate Stands			4 sets Binocular

The XSZ Series of MultiViewing biological microscopes represent great value for money. There are 3 models to choose from, each includes a complete set of good quality achromatic objectives, Abbe condenser and mechanical stage. A powerful 100W illumination source is included in the 510 model to ensure optimum viewing for multiple, remote observers.

The viewing ports give crisp, bright images of the scene under observation by the instructor. In the 510 model a useful illuminated teaching indicator is included as standard.

This microscope system is perfect for all teaching and consulting applications.

Supplied to you by:



www.gxmicroscopes.com

Due to our policy of continuous development specifications may change without notice

June 04